# Centrifuge 5415 D and Centrifuge 5415 R

Instruction Manual · Mode d'emploi succinct  
Instrucción sucinta

**eppendorf**

**Centrifuge 5415 D / R**

Instruction Manual . . . . . 5

Mode d'emploi succinct 5415 D / R . . . . . 18


Instrucción sucinta 5415 D / R . . . . . 23

No part of this publication may be reproduced without the prior permission of the copyright owner.  
Copyright© 2001 Eppendorf AG, Hamburg

# Centrifuge 5415 D


Figur 1

Figur 1

# Centrifuge 5415 R

Figur 2

# Contents

<b>1</b>	<b>Introduction</b> .....	<b>6</b>
<b>2</b>	<b>Safety precautions and applicational limitations.</b> .....	<b>7</b>
<b>3</b>	<b>Operation</b> .....	<b>9</b>
3.1	Control elements .....	9
3.2	Loading the rotor .....	9
3.3	Centrifugation with time setting .....	9
3.4	Zentrifugationszeit "AT SET RPM" (for 5415 R only) .....	10
3.5	Short-spin centrifugation .....	10
3.6	Continuous centrifugation .....	10
3.7	Switching to the rcf display .....	10
3.8	Cooling (for 5415 R only) .....	11
3.9	Opening the centrifuge in the case of a power failure .....	11
3.10	Device fuses .....	11
<b>4</b>	<b>Cleaning and maintenance</b> .....	<b>12</b>
<b>5</b>	<b>Troubleshooting</b> .....	<b>13</b>
<b>6</b>	<b>Technical data</b> .....	<b>15</b>
<b>7</b>	<b>Ordering information</b> .....	<b>17</b>

# 1 Introduction

The Centrifuge 5415 D is a non-refrigerated bench-top centrifuge and the Centrifuge 5415 R is a refrigerated bench-top centrifuge. Both are used in research and routine laboratories in the biosciences, industry, clinics and chemistry.

With their space-saving design and user-friendly operating concept, these centrifuges are ideal for everyday routine tasks.

In both devices 24 or 36 Eppendorf micro test tubes can be centrifuged simultaneously in a 45° fixed-angle rotor at 16,100 rcf (13,200 rpm).

The basic functions can be selected using the "SHORT" key and the "TIME" and "SPEED" dial. The values appear in the display in digital form.

Before starting up Centrifuges 5415 D / R for the first time, please read the instruction manual.

This sign is found on your centrifuge and on several pages in the operating manual. Texts labeled with this sign contain safety notes. Read these safety precautions before using the centrifuge for the first time.

## 1.1 Delivery package

- 1 Centrifuge 5415 D with fan cooling or 1 Refrigerated Centrifuge 5415 R
- 1 Main power cable
- 1 Instruction manual
- 1 Rotor key
- 1 Set of fuses

## 1.2 Setting up the device

Place the centrifuge onto a stable, horizontal work surface. Ensure that the working environment is well ventilated! According to the requirements of the EN 61010-2-020 standard, a safety distance of 30 cm must be observed around the centrifuge during operation. No objects whose destruction may cause damage should be positioned in this space.

Please remove the transport safety device from the rotor chamber and keep it in a safe place in case the centrifuge should be transported in the future.

**For 5415 R only:** To avoid damage to the compressor caused by incorrect transportation, wait for four hours after setting up the device before switching on.

Before plugging in the centrifuge, compare your power supply with the electrical requirements listed on the identification plate.

Fuses in device:	230 V	2.5 A, slow-acting
	120 V	5 A, slow-acting (UL fuse) (only 5415 D)
	120 V	6.25 A, slow-acting (UL fuse) (only 5415 R)

Do not operate centrifuges which have not been correctly installed or repaired.

Insert a rotor and fasten using the rotor key. Plug in the centrifuge and switch it on using the switch in the right-hand side.

The centrifuge is ready to operate when the figures appear in the display.

Please do not move the centrifuge while the rotor is spinning.

## 2 Safety precautions and applicational limitations

For your own personal safety, please observe the following safety precautions:

The rotor and the rotor lid must always be securely fastened.

The rotor must be loaded symmetrically. Opposing tubes should be of the same type and should be filled equally.

Rotors are high-grade components which are subject to extreme mechanical strain. Aluminium rotors are protected against corrosion caused by commonly-used laboratory chemicals by means of an electrolytic coating, although this protection cannot be fully guaranteed.

Aggressive chemicals, such as concentrated and mild alkalis, concentrated acids, solutions containing mercury ions, copper ions and other heavy-metal ions, chlorinated hydrocarbons and concentrated saline solutions can damage the rotors.

In the event of contamination caused by aggressive agents, the rotor must be cleaned immediately. This is particularly important to the bores.

Please clean your rotor regularly using a **neutral** cleaning agent (e.g. Extran® neutral or RBS neutral). This will protect the rotor and prolong its service life.

Please ensure that the rotor is protected from corrosion and mechanical damage. Even slight scratches and cracks can cause severe damage to the rotor materials.

**Damaged rotors must not be used in the centrifuge!**

**Before starting the device, ensure that the rotor is fastened correctly.**

Repairs must only be performed by an Eppendorf authorized service technician.

Only use original rotors and spare parts from Eppendorf.

The Centrifuge 5415 D / R may only be used for the applications mentioned in this manual.

The centrifuge must not be operated in a hazardous or flammable environment and must not be used to centrifuge explosive or highly reactive substances. The centrifuge has to be used by trained professionals.

When handling toxic or radioactive liquids or pathogenic bacteria belonging to Risk Group II (see World Health Organisation: "Laboratory Biosafety Manual"), observe national regulations. Extensive instructions regarding the handling bacteria or biological material from this risk group or higher risk groups are laid down in national and international regulations.

If such liquids are spilled in the rotor or the rotor chamber, the centrifuge must be cleaned thoroughly.

Before using cleaning and disinfecting methods other than those specified, please check with the manufacturer that these methods do not damage the centrifuge.

A liquid density of 1.2 g/ml must not be exceeded at the maximum rotational speed.

Sample tubes heat up slightly in centrifuges with fan cooling. Observe the limiting data regarding the load-bearing capacity of the centrifuge, as specified by the manufacturer.

The use of organic solvents (e.g. chloroform, phenol) may have an adverse effect on the stability of plastic test tubes.

## 2 Safety precautions and applicational limitations

2

Before being transferred from a **cold room** to a normal laboratory, the Centrifuge 5415 D / R must be run for 30 minutes in the cold room until it is warm. Alternatively, allow it to warm up for at least three hours in the laboratory, but **do not plug in the centrifuge** in order to prevent damage caused by condensation.

Lids which show visible signs of damage must be replaced.

When handling the rotor lid, please observe the resistance of the materials of construction (polypropylene / polycarbonate), e.g. the limited chemical resistance of polycarbonate to phenol.

Watch your fingers when closing the centrifuge lid.

## 3 Operation

### 3.1 Control elements

See the front view of the device, found on the fold-out page at the front of this manual:

- **TIME** – Turn time-selection dial
  - **START/STOP** – Press time-selection dial
  - **TEMP** – Temperature-selection keys (only 5415 R)
  - **SHORT** – Short-spin key
  - **FAST COOL** – Preliminary cooling run key (only 5415 R)
  - **SPEED** – Turn speed-selection dial
  - **RPM/RCF** – Press speed-selection dial
  - **OPEN** – To open the lid
  - Rotor with 24 or 36 bores for microcentrifuge tubes
  - Rotor lid with central fastening button
  - Mains power switch (right-hand side)
  - Mains power socket with fuse element (5415 D rear, 5415 R right-hand side)
  - Identification plate with power supply values (rear)
- E** Emergency lid release, for 5415 D on the base of device, with cord; for 5415 R with rotor key from right.

### 3.2 Loading the rotor

The micro test tubes with a nominal volume of up to 2.0 ml are loaded into the bores of the rotor. The test tubes are positioned in pairs opposite each other and should each contain roughly the same amount of sample.

The specifications for the maximum permitted weight per bore apply to the adapter + test tube + contents.

In order to minimize air friction and noise, the rotor should not be operated without the rotor lid.

### 3.3 Centrifugation with time setting

Switch on the power supply. The nominal values from the last run appear in the display. Load the rotor **symmetrically**. Attach the rotor lid and close the lid of the centrifuge.

- |              |   |
|--------------|---|
| <b>TIME</b>  | To modify the run time.   |
| <b>SPEED</b> | To modify the speed in increments of 200 rpm or in rcf.           |
| <b>START</b> | To start the run. The ■ symbol flashes when the rotor is running. |
| <b>STOP</b>  | To abort the run.   |
| <b>Temp</b>  | To set temperature preselection (only 5415 R)                     |

The time setting and the rotational speed can be modified during the run. The new nominal rotational speed starts two seconds after the modification has taken place.

## 3 Operation

3

The remaining time appears in the display during the run.

When the time is changed during a run, the time which has already elapsed is taken into account.

The last minute is counted down in seconds.

After the set time has expired, the centrifuge stops automatically. If the time is modified during the braking phase, the centrifuge accelerates to the nominal rotational speed.

### 3.4 Centrifugation time setting "AT SET RPM" (for 5415 R only)

When operating in the "AT SET RPM" mode, the centrifugation time counts down after reaching the preset speed.

**START** If the centrifuge lid is open and you press it for longer than 2 s, the centrifuge will switch the operating mode to "AT SET RPM" and "AT SET RPM" status will be displayed. If you wish to exit the START mode when the centrifuge lid is open, keep it pressed until the status message has been canceled.

### 3.5 Short-spin centrifugation

The SHORT key activates a short-spin run at a preset rotational speed or at maximum speed for as long as the key is held down.

**SHORT**

- If the centrifuge lid is closed and you press the key for longer than 2 s, the centrifuge will switch the mode to short spin at preselected or maximum speed.
- If the centrifuge lid is open and you press the key for longer than 2 s, the centrifuge will switch the mode from "maximum speed" to "preselected speed" and vice-versa.
- If the centrifuge lid is open and you press the key briefly, the mode you have selected will be displayed  
13 t      Short-spin centrifugation at maximum speed  
1 – 13 t   Short-spin centrifugation at preselected speed.
- a short-spin run which has already been interrupted can be continued by pressing this key, if the centrifuge has not yet come to a standstill. This may be done three times in total. After this, the lid must be opened.

### 3.6 Continuous centrifugation

**TIME** Using the keys, set the dial to above 99 min or below 0.5 min to display 00 (= ∞). The time is counted upwards in minutes.

**STOP** To end continuous operation.

### 3.7 Switching to the rcf display

**SPEED** If the Speed key is pressed, the rpm and rcf values appear alternately in the display.

## 3 Operation

### 3.8 Cooling (only 5415 R)

Using the **TEMP** keys, it is possible to set the nominal temperature between 0 °C and + 40 °C. The nominal temperature can also be changed during the run.

When the nominal temperature that has been set has been reached, a deviation in excess of  $\pm 3$  °C will be indicated by a flashing temperature display.

If the temperature deviates by more than 5 °C after the nominal temperature has been reached, the centrifuge switches off and "err 18" appears in the display.

**Fast Cool** starts a temperature-control run at a rotor-specific speed. It enables the new, preselected temperature to be attained rapidly in the rotor. (Cooling down from room temperature to 4 °C takes max. 16 minutes.) The run can be stopped automatically or by pressing the **Stop** key. An intermittent acoustic signal is emitted.

**Standby cooling:** With the lid closed, the rotor chamber is cooled down to the preselected nominal temperature before or after the run, on condition that this temperature is below the ambient temperature. Likewise, Centrifuge 5415 R maintains the temperature after a run when the lid remains closed, temperature adaptation is slow. The nominal temperature appears in the display.

If the centrifuge is not used for longer than six hours, or if the lid is not opened for this period, the refrigeration function switches off in order to protect the device. The actual temperature in the rotor chamber appears in the display. The desired temperature can be rapidly reattained via **Fast Cool**.

Remove condensation water or ice (by thawing) from the rotor chamber at regular intervals using a soft, absorbent cloth.

**Note:** For problem-free operation at room temperatures below 18 °C, a warm-up time of approximately one hour is necessary.

### 3.9 Opening the centrifuge in the case of a power failure

Unplug the device. Wait until the rotor comes to a standstill.  
(This may take up to 3 minutes!)

The emergency lid release of the 5415 D is activated by pulling a cord on the base of the device:

A small white plastic button is located in the ground plate behind the front, right-hand suction foot. (see "E" on the fold-out page at the front of this manual). Pull out this button and pull the cord vertically downwards.

#### **Emergency lid release 5415 R:**

Remove the round stopper on the right-hand side of the housing. Insert the rotor key provided through the bore that has been vacated. The centrifuge lid is released by rotating to the left.

### 3.10 Device fuses

The fuse holder of the 5415 D is located below the mains plug, the fuse holder of the 5415 R is located on the left of the main switch. These holders may be pulled out backwards and the two fuses may be replaced (see Ordering information).

## 4 Cleaning and maintenance

4

The outside of the centrifuge and the inside of the rotor chamber should be cleaned once a month using a warm detergent. Before doing so, disconnect the main power plug. Remove the rotor and clean separately. The rotor chamber should only be cleaned using a moist cloth. Only **neutral** agents may be used for cleaning and disinfecting (e.g. Extran® **neutral**, RBS neutral or 70 % ethanol/water mixture). Please ensure that no liquid enters the slit on the outlet of the motor shaft.

The rotor should be cleaned once a month using warm detergent. If the rotor can't be removed, press down one side of the rotor edge.

Check the rotor and rotor bores for residue, corrosion and mechanical damage.

The rotor is then reinserted and tightened using the rotor key provided.

All rotors are autoclavable (121 °C, 20 min.). Please check the rotor and rotor bores for corrosion and mechanical damage on a regular basis.

**5415 R:** If condensation water forms in the rotor chamber, dry the chamber using a soft, absorbent cloth.

Switch off the centrifuge after use and leave the lid open.

## 5 Troubleshooting

<b>Error</b>	<b>Cause</b>	<b>Solution</b>
No display.	No main power connection.	Plug in mains cable on both sides.
	Power failure.	Check main power fuse of device and of laboratory.
Lid cannot be opened.	Power failure.	Emergency lid release (see Sec. 3.9).
	Rotor is still running.	Wait for the rotor to come to a standstill.
Centrifuge shakes during acceleration.	Rotor not loaded symmetrically.	Stop centrifuge run and load rotor symmetrically.
	Rotor lid does not open correctly.	Stop centrifuge. Close the lid.
Centrifuge does not start.	Lid not closed.	Lift the lid and close it again.
The centrifuge does not restart in the Short-spin mode.	Safety function to prevent centrifuge from overheating.	Do not release key during short-spin centrifugation.
Flashing temperature display (only 5415 R).	Temperature deviation of > 3 °C from nominal value.	
Display <b>LID</b>	Lid latch not disengaged.	Close lid correctly, press Start / Stop or activate emergency lid release (see Sec. 3.8).
	Lid latch not engaged.	Lift up lid slightly and close it again.
<b>INT</b>	Power supply interrupted during run.	Check main plug. Repeat run.
<b>OFF</b>	Mains supply is switched off.	
<b>Error 3</b>	Errors in rotational speed measurement system.	Leave centrifuge switched on and wait for 5 minutes. Repeat run.
<b>Error 5</b>	Lid switch.	Wait for the centrifuge to come to a standstill and open the lid. Then close the lid and restart the centrifuge.

## 5 Troubleshooting

Error	Cause	Solution
<b>Error 6</b>	Drive error.	Repeat run.
<b>Error 7</b>	Overspeed.	Check whether rotor is positioned correctly.
<b>Error 8</b>	Rotor loose. Drive error.	Tighten rotor.
<b>Error 9 – 25</b>	Electronics error.	Repeat run.
<b>Error 18</b> (only 5415 R)	Temperature error. ( $\Delta \geq 5 \text{ }^\circ\text{C}$ ).	Switch off the centrifuge.
<b>Error 19</b> (only 5415 R)	Cooling circulation is overheated.	Make sure that the ventilation slit is not obstructed.
<b>Error 24</b> (only 5415 R)	Error with cooling aggregate (e.g. overheating).	After long period of operation: Allow to cool.  At the beginning of operation: Contact Service.

If the solutions suggested here prove to be unsuccessful, please contact Brinkmann Service at 800-645-3050 in USA and 800-263-8715 in Canada.

## 6 Technical data

### Centrifuge 5415 D

Power supply:	230 V / 50 – 60 Hz 120 V / 50 – 60 Hz see identification plate
Power requirement:	180 W
Max. rotational speed:	13,200 rpm
Max. centrifugal force:	16,110 rcf
Max. load:	24 x 2.0 ml Safe-Lock Micro Test Tubes or 36 x 0.5 ml tubes
Max. kinetic energy:	3,100 Nm
Max. density of material to be centrifuged:	1.2 g/ml
Ambient temperature:	2 °C – 35 °C
Max. relative humidity:	75 %
Dimensions:	Height: 230 mm Depth: 310 mm Width: 230 mm
Weight (without rotor):	8.5 kg
Acceleration time (230 V):	< 13 sec.
Deceleration time:	< 11 sec.
Acceleration time (100 / 120 V):	< 15 sec.
Deceleration time:	< 13 sec.

### Centrifuge 5415 R

Power supply:	230 V / 50 – 60 Hz 120 V / 50 – 60 Hz see identification plate
Power requirement:	300 W
Max. rotational speed:	13,200 rpm
Max. centrifugal force:	16,100 rcf
Max. load:	24 x 2.0 ml Safe-Lock Micro Test Tubes or 36 x 0.5 ml tubes
Max. kinetic energy:	3,100 Nm
Max. density of material to be centrifuged:	1.2 g/ml

## 6 Technical data

6

Technical data

Ambient temperature:	15 °C – 35 °C
Max. relative humidity:	75 %
Dimensions:	Height: 249 mm Depth: 450 mm Width: 290 mm
Weight (without rotor):	21.2 kg
Degree of pollution	2
Overvoltage category	II
Acceleration time (230 V):	< 13 sec.
Deceleration time:	< 11 sec.
Acceleration time (100 / 120 V):	< 15 sec.
Deceleration time:	< 13 sec.

Technical specifications subject to change!

## 7 Ordering information

<b>Centrifuge 5415 R</b> 120V / 50 – 60 Hz, without rotor	22 62 140-8
<b>Centrifuge 5415 D</b> 120 V / 50 – 60 Hz, without rotor	22 62 121-1
<b>Standard rotor F-45-24-11</b> aluminum, 24 bores, max. $\varnothing$ 11 mm, 45° angle, for 1.5 / 2 ml microcentrifuge tubes, with polypropylene lid	22 63 650-2
<b>Rotor lid</b> aluminium, aerosol-tight, for F-45-24-11, incl. 1 replacement sealing ring	22 63 655-3
Adapter for 0.5 ml microcentrifuge tubes and 0.6 ml Microtainers, for F-45-24-11, set of 6	22 63 622-7
Adapter for 0.4 ml microcentrifuge tubes, for F-45-24-11, set of 6	22 63 624-3
Adapter for 0.2 ml PCR test tubes, for F-45-24-11, set of 6	22 63 626-0
<b>Special rotor F-45-36-8</b> aluminum, 36 bores, max. $\varnothing$ 8 mm, 45° angle, for 0.5 ml microcentrifuge tubes, with polypropylene lid	22 63 650-2
Adapter for 0.2 ml PCR test tubes, for special rotor F-45-36-8, set of 6	22 63 660-0
<b>Accessories</b>	
Replacement lid for rotor F-45-24-11 (polypropylene)	22 66 728-9
Replacement lid for rotor F-45-36-8 (polypropylene)	22 66 729-7
Sealing rings for aerosol-tight rotor, set of 2 x 5	22 66 731-9
Rotor key	22 66 718-1
Fuse 230 V (2.5 A slow)	22 66 818-8
Fuse 120 V (5 A slow UL)	22 66 822-6

### Important:

Please use the original accessories recommended by Eppendorf. Using spare parts or disposables which we have not recommended can reduce the precision, accuracy and life of the centrifuges. We do not honor any warranty or accept any responsibility for damage resulting from such action.

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eppendorf® is a registered trademark.